

DCI PharmAssist[®] Lift Assist Solutions

Standard Horizontal Mount Lift Assist

- Linear open assist for any size manway cover.
- Compact, sealed and cleanable design.
- All stainless steel design including Belleville springs.
- Bushings and gasketing of FDA approved material.
- Standard finish is 32R_a, other finishes available upon request.
- Adjustable design.



Vertical Mount Open/Close Assist

Same features as standard assist, plus...

- Linear open/close assist specially designed for vessels with flanged top head.
- Belleville springs loaded in both directions for ease of opening and damping any overhung load.
- Compact right-angle design allows the assist to remain close to the vessel side, eliminating interference and safety concerns.



DCI, Inc. – Corporate Headquarters
600 North 54 Avenue (56303)
P. O. Box 1227 (56302-1227)
St. Cloud, MN
Phone: 320-252-8200
Fax: 320-252-0866
Web: www.dciinc.com
Email: sales@dciinc.com