



**PharMix**  
Mixers & Agitators



*guarded cage turbine type*

*propeller type*

*radial turbine type*

*axial turbine type*

*squirrel cage*

*Rushton turbine*

*dispenser type*

*bottom & side sweep*

*bottom & side scraper*



# PharMix

dcil Mixers & Agitators



- Factory operational test runs
- Pre-validation services
- Sanitary design
- Dry running, purged & double mechanical shaft seals
- Integrated controls
- Built to current GMP & ASME BPE standards
- Capabilities from 1/4 HP to 300 HP+
- Guaranteed cleaning & performance
- Full range of designs for any size, viscosity & agitation duty or service conditions

DCI, Inc.—Corporate Offices  
600 North 54 Avenue (56303)  
P.O. Box 1227 (56302-1227)  
St. Cloud, MN

Toll Free 800-276-7599  
Phone 320-252-8200  
Fax 320-252-0866  
[www.dciinc.com](http://www.dciinc.com) · [sales@dciinc.com](mailto:sales@dciinc.com)



An Employee  
Owned Company