

# Cosmetic



We serve the cosmetics and personal care industries with bulk storage vessels and equipment for mixing, blending, and batching applications. As needs change, we continue to adapt and develop custom solutions for stainless steel vessels and process tanks.

## Products

- Storage Tanks and Vessels
- Media Prep Tanks and Vessels
- Mixing Tanks and Vessels
- Blending and Batching Tanks and Vessels
- Portable Tanks and Vessels
- Agitation (PharMix)

Not only does DCI Inc. manufacture the above equipment for the cosmetic industry *but we build custom engineered equipment* that fits your company needs. Send DCI, Inc. your existing prints and specifications to get started. If you don't have existing prints please reach out with the information you do have and allow DCI Inc. to help engineer your new equipment.

## Capabilities & Customization

DCI, Inc. will customize any piece of equipment to your exact specifications. Get the most out of your equipment with our value added solutions and services:

- Material/Weld Finishes
- PharmAssist Lift Assist Solutions
- Heat-transfer Surfaces
- Passivation
- Riboflavin Testing
- PharMix Agitators with customized CIP/SIP capable seals
- Control/Automation Packages
- Integrated Systems
- Positive Material Identification
- Electropolishing
- Variety of Heat transfer Applications (channel, dimple, half-pipe, etc.)
- Catwalks, Ladders, and Platforms

## Additional Services

- Field Fabrication At Your Site
- Field Service - Repairs and Modification



DCI, Inc.  
600 North 54th Ave  
St. Cloud, MN 56303  
Phone: 320-252-8200  
[www.dciinc.com](http://www.dciinc.com)

Stay up to date with everything DCI, Inc.!

